

ACF-100 Automatic Capsule Filling Machine



Machine identity:

- Tamping filling, get better accuracy and more average weight.
- Vacuum pump using German made dry pump, oil free to gain the usage of the pump, it also easy to clean and maintaining.
- Filling weight adjustable.
- Additional large filter to gain the filter effect making the longer usage of pump and little maintaining, cleaning.
- Fills the powders that running un-smoothly, such as Ganoderma Lucidum, pollen, and royal jelly etc.
- Good quality, easy maintaining and operation.
- IQ/QQ documents are available.

Described:

- 1. Empty capsules are fed from a storage hopper to the rectifier inserted with the aid of vacuum into one of the 10 sets of carrier segment, each get 2 rows holds 18 empty capsules.
- 2. Second row for the empty capsule feeding.
- 3. Upper carrier segment goes up and pushing in.
- 4. Powder slug is inserted into of the capsule body; when there is no capsule body held at the lower holder, the powder slug is further pushed downward and collected in a receiver trav:
- 5. This station reserved for attachments to fill pellets or other granulated products.
- 6. Ejection of un-separated capsules or the up-side-down capsule caps is performed with a faulty ejector.
- 7. Upper carrier segment pushing up and then moves down.
- 8. At this station, upper and lower holders are closely matched and 10 stroking rods push the bodies into the caps.

- 9. Now the closed capsules arrive at this station are again pushed upward by the other 18 steel rods, and leave the carrier segment.
- 10. Carrier segment is vacuum-cleaner and proceeds to the capsule insertion station.

Technical Data:

Model No	ACF-100
Capacity (caps/hr)	58,000~62,700
Capsule Size	#00, #0, #1~#4
Power Supply	220/380V 3Phase
	50/60Hz
Power Consumption	20.6 A
Main drive	2 HP
Auger drice	1/2 HP
Vacuum Pump	2 HP
De-duster	3HP*2 Pcs
Machine Weight	1200 KGs
Dimension (L×W×H)	983×1201×2013 mm