Optional item



The Pellet device is for filling pellets into hard gelatine capsules. It is designed to minimize weight deviation so the supply of pellets can be easily adjusted to achieve the target weight.



This provides the powder automatically to the machine's hopper

Both Auto & Manual mode are available so the machine is ideally equipped for an unmanned production environment.

Pellet filling device **Change part**

The machine can be configured for all capsule sizes with a range of change parts available.

The Capsule Filler is designed for ease of changeover by competent operators.



Other option

O HMI system IQ, OQ

Dust collector Calibration

O Capsule polisher & sorter

Powder loader

Capsule loader

The Capsule Loader is designed to elevate the capsules to prevent capsule separation and crushed capsules during

Automatic operation is controlled by a sensor.





The above specification may be subject to change without notice.

O Why choose CAPLER™ capsule filler?

- **Good Performance**
- Ergonomic Design
- **Technological Innovation**
- √ User Friendly Design
- Rapid Change Over & Cleaning


Pharm Tech Korea

Pharmaceutical Machine

Distributor / Representative

#21 Majun-Dong, Seo-Gu Incheon City, Korea. 404-820 +82-32-564-6851~4 +82-32-564-6855

Sales@pharmtechkorea.com Web site http://www.PharmTechKorea.com http://www.PharmTechKorea.co.kr



HMI system


The operation screen of the **CAPLER**™ is a fine example of 'User-friendly 'design.

The software ensures the optimum capsule quality with its production data management and automation control. The product data and settings can be stored and re-called to aid uniform production.

Operation

Easy and organized screen composition will bring satisfaction even to the beginner.





All the production data can easily be saved, providing data for product traceability and machine history. Therefore retrospective information about when, who, what, how and why can be investigated and reported.

Report

Group management

A management group can be created, deleted and modified by the authorized manager to control user access levels.

Group List		c	Create Group		nfo Delete Gr	Delete Group	
No	Group Name		Descr	ription	Create Date	1	
1	4000		ant	100	07-14-2007 12:00		
2	Guest		Lev	re(0			
3	Supervisor teset level	2 Ower Gross			10-06-2005 18:30		
	40.000.000.00	Group Notes			19.42.0994.2028		
		Description 1					
		Assess Printege Check Lint					
			Access Phiologe Tree	Chert Cade 40			
			ar All Authority				
			Operation	******			
Committee of the Commit		ar Product 808 608 608				- 4	
Access Privilege List		- Safety en mo me					
No	Access Privilege		Report Management	212 000 000			
1	Time Reset	+c System ************************************			Reset Button		
	71009,750,441		Product Print	444 500 500	CAST COMMON		
			Print Screen	\$100 may \$100.			
			QUIT	100 000 000			
			Drive On	114 mm mm			
		400					
				i i			
			Circum	Onne			
						iel	



In the event of a machine stop by an alarm, it is easy to check the related alarm at a glance. The diagnostic display gives reason and solution for the activated alarm so the operator can quickly resolve the issue without engineer intervention.

> Safety message



Powder tamping

The powder will be compressed with 5 steps of tamping, the compressed powder will be inserted into the capsule at 6 step. PTK's new design of scraper will reduce the weight deviation.





Closing & Discharge

The Capsule closing and Discharge assemblies are designed for easy detachment and fitting.

This makes Maintenance and Cleaning very convenient for the user. The discharge parts are designed to prevent a vortex during capsule discharge.

Filling combination



+ Powder



+ Powder

+ Pellet



+ Powder



Pellet

+ Tablet





Tablet + Tablet + Tablet

Turn table

Precise design is applied to the Turn table, (segments and segment holders). The turn table is rotated by Index movement. The precision of segment is improved by the usage of a straight





Capsule feeding

The empty capsules are provided to the segment through the feeder. The feeder is designed for smooth feeding with a race way and side pusher. The simple construction of the race way and side pusher makes them quick and easy to change.

Fully compliant with requirements of FDA21 CFR Part11

Automatic Capsule FIlling Machine_



Pharm Tech Korea 's PF range of Capsule filling Machine embraces the concept of simplicity by good design. With logical operation and high quality components we will not fail to impress. Contact us to experience the innovative **CAPLER**™ concept!





Pharm Tech Korea

Best supporter to you

