



SHH-QCL60 Automatic Vertical Ultrasonic Glass Ampoule Cleaning Machine



SHH-QCL series automatic washing machine consist of infeeding section, ultrasonic treatment section and final washing section, the containers are going through those sections for completion of the whole washing processing by Ultrasonic treatment, prewashing and final washing steps, the cleaning quality of the container washed meets any Pharmacopoeia and being able to validate.

Working principle:

Bottles are input through conveying belt, slide into the pool after infusing water, then are washed through ultrasonic about one minute. The screw separates bottles and sends bottles to lifting device. The screw is designed according to best curve, make bottles conveyed from screw to lifting device smoothly. Lifting device lifts bottles smoothly from water and conveys to manipulators of rotary drum. Manipulators clamp bottles and turn 180°, in order to enter following washing station. First step, high-pressure circulating water flush the exterior wall of bottles, then the first group of needles insert into bottle and flush circulating water, the second group of needles flush circulating water, the third group of needles flush compressed air, discharge residual circulating water. The fourth group of needles flush inject water, the fifth, sixth group of needles flush compressed air, discharge residual water at most in order to dry.

Finally flush compressed air to exterior wall of bottles. Manipulator turnover bottles and sent to output device, then conveyed to tunnel or output plate through output star wheel. The whole washing process is finished

Features:

1. Wide range of different vial sizes from 2-50 ml, conventional format change
2. Ultrasonic treatment
3. High efficiency Washing procedures
4. Optional automatic hood lifting system
5. Glass tissue discharge device
6. Ultrasonic power level online display, recording, saving, printing functions
7. Water tank heating system
8. Low pressure alert
9. Use of the FDA guideline listed material

Technical parameters:

Model	SHH-QCL40	SHH-QCL60	SHH-QCL80	SHH-QCL100	SHH-QCL120
2-5ml	9000	18000	24000	30000	36000
7-10ml	9000	18000	24000	30000	36000
15ml	7200	14400	21000	21000	27000
20ml	6600	13200	19000	12000	27000
30ml	6000	12000	12000	9600	N/A
50ml	4200	N/A	N/A	N/A	N/A
100ml	N/A	N/A	N/A	N/A	N/A
Power supply	380V/50Hz/3Phases				
Power	15KW	15KW	15.5KW	15.5KW	15.5KW
Compressed air	≤0.01%	≤0.01%	≤0.01%	≤0.02%	≤0.02%
Compressed air consumption	0.3-0.5Mpa				
WFI temperature	≥35 m3/h	≥40 m3/h	≥60 m3/h	≥60 m3/h	≥60 m3/h
WFI consumption	40-60℃	40-60℃	40-60℃	40-60℃	40-60℃
PW pressure	0.30-0.5 m3/h	0.35-0.65 m3/h	0.40-0.8m3/h	0.40-0.9m3/h	0.50-1.0 m3/h
Dimension	1795*1894	2119*2306	2230*2381	2230*2381	2230*2381
	*1300	*1300	*1300	*1300	*1300
Feeding height	1000mm	1000mm	1000mm	1000mm	1000mm
Weight	2100KG	2700KG	3100KG	3110KG	3110KG