




Tablet Dispenser PP-03



The counter PP-03 is a reliable, high performance machine, designed in conformance with current and proposed GMPs. Having a small footprint, they are effective in maintaining the integrity of a clean room environment. Few change parts are needed for different container sizes (only the nozzle). All product contact parts are made of type 316 stainless steel. And all contact parts can be easily disassembled for cleaning. It is applicable to horizontal and vertical pouch packaging line.

Products (A) are loaded into the hopper (B) manually or automatically (option), and then fed onto the 1 set of three-grade vibratory plates (C). By adjusting the vibration frequency of the vibratory plates, the products can be allowed to move through 8 channels of the plates continuously and in single file. At the end of channels the products fall into 8 detection tunnels (D), where the products are counted with the help of combination of sensors (E) and control system. The counted products drop through the funnel (G) and the nozzle (I), then into the bottle. When the preset count is reached the tunnel gate (F) closes. After the bottle is filled, the bottle positioning cylinders work to release the filled bottle and get the next empty bottle positioned under the nozzle.



Technical parameters:

Item	Specification	Remark
Max. output	3,000 pcs/min	Dependant on product and container size
Fill qty.	1~1000	
Product size	Thickness:3~10mm Dia:3~13mm Length:3~21mm	Suitable for tablets, capsules, pills, caplets, soft gels etc.
Container size	Dia.:20~100mm Height:40~240mm	
Power supply	AC220V 50/60HZ	
Wattage	2KW	
Air pressure	0.5~0.7Mpa	
Air consumption	80L/min, clean air	
Measurement(L×W ×H)	Approx. 1717mm×806mm×1757mm	
Weight	Approx. 600kg	