



## *XHL-DXD-series Pharmaceutical Dry Powder Filling Machine*



This machine absorbs the advanced experience of the bag-filling packaging machine at home and abroad, introduces the latest Western technology, adopts PLC programmable control, combines the advanced touch screen man-machine interface system, and is an ideal metering equipment of the large-scale bag-filling packaging machine at present. The machine is reasonable in design and novel in appearance. Use Standard Strip cover, change the packing, the machine can be used powder and water, according to the requirements, can be used drag-type three-side seal and four-side seal.

### **Packing process**

1. Film shelf
2. Shaping
3. Packing film jig
4. Bottom sealing
5. Checking color base
6. Vertical sealing
7. Filling
8. Top sealing
9. Shearing
10. Towing bag
11. Packing Product

## Technical parameters

Product Model	DXD-110	DXD-180A	DXD-180B
Applied to material	Powder, liquid, viscosity		
Filling heads	1	1	2
Packing speed	60-80bpm	40-80bpm	80-120bpm
Suit to bag width	50-110mm	80-180mm	50-90mm
Suit to bag height	50-130mm	100-220mm	50-110mm
Packing volume range	2.5-60ml	2.5-380ml	2.5-190ml
Roll film out-diameter	≤500ml		
Power consumption	1.9kw	5kw	5.5kw
Measure error	≤±2%	≤±1%	≤±1.5%
Net Weight	700kg	950kg	1200kg
Overall dimension(mm)	2060*620*1600	3260*1020*2390	3260*1020*2390