

Dry Powder Vial Filling Semi Plug Pressing Machine XQBZ-30



Application:

1. It is widely used in precision, quantitative loading / dilution of pharmaceutical, pharmaceutical, food, chemical, cosmetic and other industries.
2. Pre-chromatography preparation, pharmaceutical QC and R&D, environmental monitoring, mineral and blood analysis, forensic identification.
3. Used for precise quantitative loading / dilution of blood, biological reagents, oral agents, injections, advanced cosmetics, flavors, fungicides and other liquids.

Characteristic

This equipment has been tested to meet the requirements of GMP certification.

1. The intelligent operating system is adopted.
2. The liquid filling part adopts high purity ceramic pump, which controls the injection volume of each lap by adjusting the angle of the pump head, and has extremely high measuring accuracy and repeated consistency.
3. 304 stainless steel appearance, desktop structure, practical and convenient.
4. Riveting part adopts special material coating hardening treatment, which has a long service life.
5. The equipment is fast installed by mold, and it can be replaced conveniently between different bottle types.
6. This device adopts the German sensor to achieve the functions of no filling without the bottle, no pressing without lids.

Technical parameters

Model	XQBZ-30				
Pump head	High purity ceramics				
Filling accuracy	$\leq\pm 0.4\%$				
Filling range	0.1-2ml	0.5-5ml	1-10ml	2-25ml	5-60ml
Manufacturing speed	25- 35 psc / minute				
Size of Host Machine	1500*1000*1000(mm)				
Power	AC 220V 50-60Hz 1500W				
Air supply	0.4-0.6Mpa				