



Automatic Inline Type Pharmaceutical 4 Heads Syrup Filling and Sealing Machine TALF-2/4



It is widely used in daily pharmaceutical syrup, chemical, food, beverages, medicine and other paste or low viscosity liquid(glass or plastic bottle with aluminum cap)

Description

Automatic syrup filling machine depends on the different products and various filling requests of the customers, for example, high viscosity or sticky cream, paste, such as, tomato paste, ketchup, jam, sauce, butter, chocolate paste, cosmetics cream, etc., low viscosity liquid products, such as, juice, beverage, shampoo, lotion, laundry detergent, dish detergent, wine, soybean sauce, vinegar, edible oil, cooking oil, drinking water, etc., and liquid chemical products, liquid pesticide, liquid fertilizer, liquid medicine, etc.,

1. With organic glass chamber
2. With feeding bottle plate-form and collection bottle plate-form
3. English operation interface
4. Machine moving by trundle
5. Voltage: 380V 50Hz
6. If material less than bottom of hopper, machine will automatically stop working
7. Machine with filling hopper

Products Parameters

Model	TALF-2	TALF-4
Filling nozzles	2	4
Voltage	110V/220V/380V, 50HZ/60HZ, 1000W	
Filling range	10-120ml, 25-250ml, 50-500ml, 100-1000ml (Optional)	
Filling speed	20-40 bottles/min	40-80 bottles/min
Filling precision	$\leq \pm 1\%$	
Air pressure	0.4-0.8MPa	
Weight	750Kg	850Kg
Dimention(L*W*H)	230*90*180cm	250*90*180cm