

## XGS30(12/20) Bottle washing machine



## Main usage:

This production line is suitable for 50ml,100ml,250ml and 500ml glass I.V solution glass bottle washing, filling, plugging, etc.

## **Performance & characteristics:**

- 1.make washing ,filling and capping in integration can greatly decrease the prophase investment and anaphase service fee.
- 2.ultrasonic cleaning can avoid causing dead angle and breakage bottleneck with fur brush washing.
- 3.water is reclaimed using at a lower rank ,its pressure and time can set according to real demand, whereby decline the consumption of water
- 4.unique water tightness structure can make the rough and fine washing zone in totally disconnected.
- 5.for high autoimmunization level and photoelectric detection it can realize "no bottle ,no filling and no capping".

## Main technical parameter:

mani tooninaa paramotori	
capacity	100-300b/min
quantity of bottle-entering per time	30(12.20)bottle
gross powder	3-9 kw,380v 50hz
specitification	50,100,250,500ml
Ion water consumption	~80 ml/bottle
distilled water consumption	~100ml/bottle
GIP/SZP	125°c~30min
noise	max 75db
compressed air	0.4-0.8mpa 4m3/h
body material	304
point of adherence with liquid medicine	3161 or special silicon tube