



XLP Clinical bottle washing machine



This machine is applied in clinical bottle and dropping bottle cleaning use to turn the drum frame, and provide a dish to provide the bottle, and be sent by stir the round bottle, turn the drum intermittent into the bottle, use again back to use the water (Or running water), go to leave the son water (or distilled water) and decontamination compressed air to hand over for flush to complete the bottle-water.

Technical Parameters:

Production capacity	3600 - 6000 bottles/hour
Applied bottle	50 - 500 ml, Glass bottle
Recycling water consumption	0.5 m ³ /h Pressure: 0.25 - 0.3 MPa
Leave the son water consumption	0.5 m ³ /h Pressure: 0.25 - 0.25 MPa
Purified compressed air consumption	0.5 m ³ /min Pressure: 0.25 - 0.25 MPa
Power	2 KW
Power supply	AC,380V/50HZ 3 phase or customization
Overall size	2000 x 1400 x 1500 mm