



ZP29 Rotary Tablet Press



This machine is a kind of double press automatic rotary, press continuous machine, press granulate material to tablet.

This machine main used for pharmaceutical tablet production, and chemistry, food, electric etc industry. It is closed completely, material is stainless steel, inner table is stainless steel, can keep surface luster, and prevent pollution, meet the GMP standard.

There is a transparent window, can watch the status of press tablet, side plate can open, easy to clean and maintain. Layout of all the controller and operation parts is reasonable.

Use frequency converter adjust speed to electric adjust speed, easy to operate, rotary stable, safety and accuracy. Equipment overload protection device, when the pressure is overload, stop machine automatically, touch button and screen display. Traditional system enclosed in the tank below the machine body, and is completely separate individual components, will not produce oil pollution, and the transmission system components infiltration in the oil sump, easy to heat and wear. Machine has dust collect device, can collect the dust in press room.

Technical Parameters:

The number of punches	29
Max. Pressure (KN)	100
Max. Diameter (mm)	22
Max. filling depth (mm)	15
Max. Thickness (mm)	8
Rotary table speed (rotary/min)	14-36
Capacity (pcs/h)	120000
Power (KW, rotary/min)	4kw
Dimensions (mm)	1500*1300*1650
Weight (KG)	1850