

SED-1-1-1XGS-G Penicillin Bottle Liquid Filling And Capping Machine



Product Overview of Filling and Capping Machine

1. It is widely used in precision, quantitative loading / dilution of pharmaceutical, pharmaceutical, food, chemical, cosmetic and other industries.
2. Pre-chromatography preparation, pharmaceutical QC and R&D, environmental monitoring, mineral and blood analysis, forensic identification.
3. Used for precise quantitative loading / dilution of blood, biological reagents, oral agents, injections, advanced cosmetics, flavors, fungicides and other liquids.

Feature

1. Automatic dual-channel bottle loading, filling, capping and discharging.
2. The machine is made of stainless steel 304.
3. This equipment adopts Japanese sensor, Taiwan precision cam for the accurate positioning;
4. The filling part adopts the precision ceramic pump developed independently, and the repeated accuracy can reach $\pm 0.3\%$.
5. Sensor recognition to realize the function of stopping filling without bottle.
6. The size of the machine can be customized according to the customer's need.

Technical parameters

Model	SED-1-1-1XGS-G
Capacity	40-45 pcs/min
Power	AC 380V 50Hz
Size	1400mm×1200mm×1350mm