



湖南汇一制药机械有限公司

HUNAN FE PHARMACEUTICAL MECHINERY CO.,LTD

## QBX *Ultrasonic Washing Machine*



### **Main Application:**

This machine can finish 50ml- 500ml glass bottles washing .

### **FUNCTION FEATURE:**

- ◆ Newly designed structure , small overall size .
- ◆ It is easy to maintain and overhaul that the driver is on the side of equipment.
- ◆ Ultrasonic rinsing can avoid the dead spot caused by the hairbrush , and the hairbrush pieces the bottleneck.
- ◆ Realizing the integration of machinery, electricity, gas. With high degree of automation . Both water pressure and the time of flushing can be set in accordance with the actual demand . Flushing bottle water can be recycled and filtered which can be used degradly. So it reduces the water consumption.
- ◆ Adopting crawler to carry bottle, it is a long time to empty water after charging distilled water. Without residual water left in bottles.

### **Main Technique Parameters:**

- ◆ Production Ability: 200-600 bottles/min
- ◆ Specification : 50ml. 100ml, 250ml, 500ml
- ◆ Power: 7-12kw
- ◆ Water consumption : common water :200ml/bottle
- ◆ Washing temperature : 40-50°C
- ◆ Noise : max75db
- ◆ Body Material: 304