

TBJ-I Type Sticking Mark Machine



Description: This machine is used to stick mark circumferentially to column-model vessel, it is a mating product of ZP-I type omni-auto-counting and bottling line of pill, tablet, soft and hard capsule of our company, it can be used stand-alone and also can be used other matching product line.

Working principle

Un-drying-glue tags are labeled in the circumferential mother strap disconnectedly in advance, tags will be peeled off automatically through photoelectrical sensor and serving generator system, and stick on the outside surface of non-column-model vessel, and thus, it finishes sticking mark process.

Technical characteristic

1. Accurate position of sticking mark, reliable, it can realize sticking mark half-circumferentially or whole-circumferentially.
2. Full function of printing data, clear writing, and position adjustable.
3. It can add up and display team output and gross per batch.
4. Simple regulation, wide adaptive range.

Technical parameter

1. Output: 50-100 bottles per minute
2. Qualification of bottle: ϕ 20-100mm (diameter) 30-150mm (height)
3. Mode of sticking mark: half circle, whole circle
4. Character of sticking mark: opacity and translucence
5. Total power: 1kw
6. Power source: three-phases AC 380V
7. Air pressure: 0.4MPa consuming gas tolerance: 0.05m³/min
8. Overall dimensions: 1800×800×1400mm (L×W×H)
9. Weight of whole machine: 150kg