

Tablet Counter Machine JF-16



Counting machine can count and pack of solid particles, such as tablets, capsules (transparent or non-transparent hard capsules, soft capsules), pills from 3 to 22mm wide into bottles cans

Introduction

Softgel Counter machine is able to combine with the vial distributing machine, cotton stuffing machine, desiccant inserter, cap screwing machine, electromagnetic induction foil-sealing machine and self-adherence labeling machine to form full-automatic high-speed bottled production line.

Liner tablet counting machine use the most advanced vibration-type multi-channel cutting, computer control, system self-check, fault indication alarm, automatic shutdown and other advanced technologies

Multi-channel tablet counting reliable smooth operation and low noise

Main Features

- 1. Strong compatibility, high accuracy, wide applicable range and adopting the dust-resistance photoelectric detection technology. It is realizing the counting and split charging for the tablet including the shaped tablet, hard capsule, soft capsule (transparent and nontransparent) and pill.
- 2. Softgel Counting Machine oscillatory type material loading, without damaging the drugs;
- 3. Patent: Bottle separation and switching device, will not damage the pills.
- 4. High intellectualism, multi detection and alarming control functions like never missing the count for vials, self-check against failure:
- 5. Multi-channel tablet counting has smart structure, stability, replacing parts infrequently, small occupation area, low cost;

- 6. Convenience in replacing the variety, automatic rise and fall, convenient and fast disassemble without any tools;
- 7. Customer required language display control panel, storing the parameters of 10-20 varieties, particular convenience in adjustment and use;
- 8. Full automatic operation, intelligent joint control between the last and next processes, perfect coordination, without person attendance, saving the labor;
- 9. Low price, particularly suitable for the manufacturers of small batch and many varieties;
- 10. Control chip is the original imported from USA, and the photoelectric sensing element is introduced from Taiwan.

Technical Parameters

Item	JF-16
Applicable target	#00-5 capsule, soft capsule, \varnothing 5.5-20 tablet, shaped, sugarcoated tablet and 3-20 pill
Production capacity	20~80 bottles/min
Description for production capacity	Capacity is subject to the type and size of the bottles and raw materials, as well as the quantity filled per bottle.
	BPM=Bottles per minute
Working range	2-9999 grains (tablets) adjustable
Power consumption	220V/50HZ/ 1.2KW
Size(mm)	1360×1650×1650mm
Filling error	≤ 2‰ bottles
Material	316L+304
Control System	MCU HMI
Pressed Air Supply	0.6 Mpa
Weight	550Kg