



B series Rotary Tablet Press



Double-pressing type and single-sided tablet discharging. It uses ZP punch to press granular raw materials into round tablets and special-shaped tablets of various specifications.

Features

1. With a device that ensures the stability of the difference in weight accuracy.
2. The stainless steel outer casing is completely closed. All the parts in contact with the medicine are made of stainless steel or treated by the surface, non-toxic and corrosion-resistant.
3. Turntable surface after special treatment, can prevent cross contamination.
4. All four sides of the tablet room for transparent plexiglass, and can be opened, easy internal cleaning and maintenance. The interior is equipped with safety lighting.
5. Variable frequency step-less speed regulating device, easy to operate, safe and reliable.
6. The hydraulic system design is reasonable, which can maintain the stability of the hydraulic system.
7. Equipped with overload protection device, when the pressure overload, can automatically stop.
8. With the function of pre-compression and main compression, which could improve the tablet quality.

Technical Specifications

Model No.	ZP35B	ZP37B	ZP39B	ZP41B
Dies(sets)	35	37	39	41
Max.pressure(kN)	80			
Max.pre-pressure(kN)	10			
Max. dia. Of tablet(mm)	13 (special-shaped 16)			
Max. depth of fill(mm)	15			
Max. thickness of tablet(mm)	6			
Turret speed (r/min)	10-36			
Max.production capacity(pcs/hour)	150000	159840	168480	177120
Motor power (kW)	3			
Overall size(mm)	1100×1050×1680			
Net weight(kg)	2300			