

KGL-W mini vial capping machine



KGL-W mini vial capping machine can seal the al cap onto the vial for pharmaceuticals, bio-chemicals, etc, after filling and coppering, the filled vial should be sealed by al cap. The mini type capping machine can be applied for sealing 2-50ml vial; it is the ideal equipment in research institute or labs.

Main technical parameter:

Vial size	2-50 ml
Productive capacity	10 - 30 vial /min.
Capping station	1
Qualified capping	≥99%
Vial breaking	≤0.1%
Speed reducer	120w
Vibrator power	40~120 w
Outside dimension	900×460×900 (mm)
Weight	120 kg