

NQF800 Deblister & Retrieving Machine



Defects, during medicine production, can cause capsules to be rejected from the capsule filling machine or Blister Machine. Pressing out the capsules from faulty blister cards and separating the capsule from the powder, without cross-contamination, can be accomplished by our machine: The Deblister and Retrieving Machine NQF800.

This is a multifunction machine that can automatically press out the capsules from rejected blister packs, separate the capsules and retrieve the powder, allowing the medicine to be recycled and reused. It is applicable for capsule No. 00#-5#.

Technical Parameters

- Max. Output: 800 capsules/min.
- Machine Power: 0.38kw
- Machine Dimension: 800×480×1500mm
- Machine Weight: 150kg