

## ZPW-15D/17D/19D Rotary Tablet Press/Pill Make Machine for Pharmaceutical



This machine is a single-press type, continuous automatic tablet press for pressing granular raw materials into tablets. It is mainly used in pharmaceutical industry and also the chemical, food, electronic, plastic and metallurgical industries.

## Operating principle:

- 1. The outside part of the machine is fully-enclosed, and it is made of stainless steel, meeting the G.M.P requirement
- 2. It has transparent windows so that press condition can be observed clearly and the Windows can be opened. Cleaning and maintenance is easier
- 3. The machine can press not only round tablets but also different geometrical shape tablets, double-layer and annular tablets, and the tablets may have the impressed letters on both sides
- 4. All the controller and devices are located in one side of the machine, so that it can be easier to operate
- 5. An overload protection unit is included in the system to avoid the damage of the punches and apparatus, when overload occurs
- 6. The machine's worm gear drive adopts fully-enclosed oil-immersed lubrication with long service-life, preventing cross pollution.

## Technique parameter:

Model	ZPW-15D	ZPW-17D	ZPW-19D
Dies (sets)	15	17	19
Max. Pressure (kn)	60	60	60
Max. dia. of tablet (mm)	22	20	13
Max. Depth of fill (mm):	15	15	15
Max. Thickness of largest table (mm)	6	6	6
Production capacity (pcs/h)	33000	40000	42000
Turret speed (r/min)	37	37	37
Motor power (kW)	2.2/960, 380V/50	2.2/960, 380V/50	2.2/960, 380V/50
Overall size (mm)	615x890x1415	615x890x1415	615x890x1415
Weight (kg):	700	700	700