

ZP-5/7/9 Rotabry Tablet Press Machine



The machine is a single-press type, continuous automatic tablet press for pressing granular raw materials into tablets. It is used for small-scale production.

Feature:

- A. The machine is compliance with the GMP requirement on the production of pharmaceuticals.
- B. The compression parts is entirely separated from the driving mechanisms
- C. All spare parts to be in contact with the pharmaceuticals are made of stainless steel or the surface being treated specially for non-toxicity & corrosion proof
- D. The compression area is enclosed with transparent organic glass. To prevent contaminates, the compression area can be specially treated.
- E. All surfaces being smooth & easy to clean.
- E. Machine speed, the depth of filling, the thickness of the tablet could be adjusted.
- F. The mechanical buffer unit will prevent punches and apparatus from overload damage.

Parameter:

Type	ZP-5B	ZP-7B	ZP-9B
Dies	5	7	9
Max.pressure(Kn)	60	60	60
Max.Dia.of Table(mm)	20	20	20
Max.Depth of fill(mm)	15	15	15
Max.Thickness of Table(mm)	6	6	6
Turret speed(r/min)	30	30	30
Production Capacity(pc/h)	9000	12600	16200
Moto(Kw)	2.2	2.2	2.2
Overall size(mm)	480*490*1100	480*490*1100	480*490*1100
Net weight(Kg)	260	260	260