



Pharmaceutical Rotary Tablet Press Machine ZPD15



The technological process of tablet pressing machine includes overfilling, quantification, prepressing, main pressing forming, tablet production and other procedures. The upper die hole of the machine is loaded with the upper punch rod, the middle die is loaded with the middle die, and the lower die hole is loaded with the lower punch rod.

Technical Parameters:

- Model No.: ZPD15
- Maximum tablet thickness: 20
- Maximum tablet diameter: 40
- Max. Tablet Diameter: 40mm
- Capacity: 13500 pcs/h
- Power: 7.5kw
- Weight: 2100kg
- Dimension(L*W*H): 1200*950*2000mm