



HM SC series Rubber Stopper Washing Machine



Model: HM SC-U10/20/30/40/50/60/80/100/120/150/180/200

Working principle: vacuum feed stopper into vessel, powerfully spray, low-speed roll and wash, silicone, steam sterilize, dry stopper, then set interlocks free to discharge.

Usage: wash, siliconization, sterilize and dry rubber stopper

Suitable: rubber stopper for vial, infusion and oral

Range: sterilization temperature 121°C; Vacuum ≤ -0.09 Mpa

Capacity: ($\phi 20 \times 8.5$ stopper) 10,000-200,000pcs/batch

Vessel: 100/250/280/400/450/600/800/950/1100/1400/1650/1800L

Cycle time: 3-5hours/batch, including 40-50mintues for sterilization