

Rotary Filling And Cover Pressing Machine XQZGD-12B



Application:

1. It is widely used in precision, quantitative loading / dilution of pharmaceutical, pharmaceutical, food, chemical, cosmetic and other industries.
2. Pre-chromatography preparation, pharmaceutical QC and R&D, environmental monitoring, mineral and blood analysis, forensic identification.
3. Used for precise quantitative loading / dilution of blood, biological reagents, oral agents, injections, advanced cosmetics, flavors, fungicides and other liquids.

Characteristic:

This equipment has been tested to meet the requirements of GMP certification.

1. Intelligent operating system with touch screen is used.
2. The liquid-dividing part adopts high-purity ceramic pump, which controls the liquid-filling quantity of each circle by adjusting the angle of the pump head, and has high metering precision and repetitive consistency.
3. 304 stainless steel appearance, desktop structure, practical and convenient.
4. The machine adopts bottle mold fast loading method, different bottle type can be changed quickly and conveniently.
5. This device adopts German electric eye detection to achieve the function of bottling without filling.
6. Pneumatic components adopt Taiwan brand.

Technical parameters

Model	XQZGD-12B				
Pump head	High purity ceramics				
Filling accuracy	$\leq \pm 0.4\%$				
Filling range	0.1-2ml	0.5-5ml	1-10ml	2-25ml	5-60ml
Manufacturing speed	0- 50PCS / minute				
Size of Host Machine	500*550*650(mm)				
Power supply	AC 220V 50-60Hz 280W				
Air supply	0.4-0.6Mpa				