



BPAD 2000II Desktop Electric Tablet Capsule Counting Machine



The machine counts and fills medicine automatically but bottles should be positioned manual. The machine adopts CCD sensitive technology and PLC controls to ensure high resolution, fast response speed and accurate counting. If the kind or filling quantity of medicine is changed, you only need to do simple setup on the touch screen. It is easy to clean and maintain due to its compact structure. Besides, it is low noise and power saving. The machine uses touch screen to realize human-machine conversation. It runs steadily and reliably and stops counting automatically once lack of medicine or no bottle under the funnel occurs. Materials in contact with medicament conform to requirements of GMP It is suitable for counting and filling a small batch of medicine in pharmaceutical, food and chemical industries and medical units.

Main technical data:

Type	Size of medicament	Size of bottles	Output	Power supply	Power	Machine dimensions (WxDxH)
BPAD2000II	TabletØ5-12mm	20-300ml	About 800-1000pc/min	AC200V 50Hz	0.35Kw	900x600x650mm
	Capsule:00#-5#					