

## Capsule Polishing Machine DXP-ABC



The capsule polishing machine is applied to the principle of brush polishing, and removes the powder attached to the surface of the capsule, make the surface of is clean and bright, and it is essential equipment for the post-production processing of capsules and tablets. The polishing process of the machine is a closed structure. The parts in contact with the drug are made of corrosion-resistant materials, are easy to disassemble and clean and meet the requirements of the pharmaceutical production management regulations (GMP) implementation guidelines for pharmaceutical equipment.

### Main Feature:

1. Small size, beautiful appearance, adjustable height, and angle.
2. It can be connected to any type of capsule filling machine.
3. The main shaft adopts fast detachable brushes and bearings. The brushes and bearings can be easily removed.
4. High-quality stainless steel is used in the polishing room.
5. The brushes of different sizes can be changed to meet the needs of different capsules.
6. There are safety devices. The speed is controlled by the frequency converter.
7. It can automatically filter out lightweight, empty or scratched capsules.

### Technical Specifications

Product Model	DXP-A	DXP-B	DXP-C
Production Capacity	3000-7000 Pcs/Min	150000 Pcs/h	150000 Pcs/h
Motor Power	250W	180W	180W
Voltage	220V	220V	220V
Machine Weight	45KG	40KG	40KG
Product Dimensions	950*400*1100MM	930*400*1150MM	930*400*1150MM