

DJL-24A Electronic Counter



The equipment counts and fills medicine automatically.

Capsule, soft gels, normal tablets, special shape tablets, coating tablets, pills and most solid granules are fed from hopper to counters by quiver to prevent material breakage, through the high anti-dust sensor technology and adjustable baffles at the outlet of loader into the bottle.

High anti-dust sensor technology ensures counting works accurately with good stability in heavy dust conditions, while baffles control the filling speed of material.

Application:

For capsule, soft gels, normal tablets, special shape tablets, coating tablets, pills and most solid granules. Suitable for filling in plastic bottle with square, round, flat-square and special shape bottle.

Technical parameters:

Bottle Sensor	To detect bottle in conveyor
Dimension (L x W x H) mm A	1800x1400x1680
Dimension (L x W x H) mm B	2000x1600x1680
Weight (kg) A	400
Weight (kg) B	460
Working Range	lebih dari atau sama dengan 2 (adjustable)
Counting system	Nantong Hengli
Speed Control	T.W.T
Capacity (Bottle per minute)	20-90
Power Supply	220 V ; 50 Hz
Power Consumption	1.5