

PHARMA MACHINERIES



SINGLE ROTARY TABLETTING MACHINE GMP MODEL



SINGLE ROTARY TABLET PRESS

DESCRIPTION

SALIENT FEATURES FOR STD & GMP MACHINE

- Flywheel and Handwheel provided.
- Friction plate is provided on cam plate to minimise wear & tear of cam track which also facilitates its easy replacement.
- Re-designed Variable speed pulley.
- Separate Centre Pillar of Meehanite grade casting.
- Turret of Meehanite grade casting.
- Very rigid Gun-metal Feedframe with sufficient height to avoid spill over of powder.
- All Columns, Nuts and Flywheel rim are Electroplated.
- Central Oil Lubrication system with handpump.
- Pressure gauge provided to check compression load.
- Excellent finish, Accessibility for ease of cleaning, operation and Maintenance.
- Oil groove on top face of turret for collection of excess oil.
- Optional equipments like Dust Extractor unit and Deburring / De-dusting unit can be attached with the machine.
- Overload pressure release mechanism to give protection against excess pressure.
- Cam tracks and pressure rolls made out of alloy steel for maximum wear resistance and long life.
- Worm is made up of Phosphor Bronze.
- RDD3-23 station machine can compress D.S. tablet upto 19mm (3/4") in length.
- Turret with Keyslots essential for shaped tooling.
- Maximum Turret RPM for STD Machine is 30.

ADDITIONAL SALIENT FEATURES FOR GMP MACHINE

- Machine is totally enclosed, i.e. Top half of the machine (Compression zone) with S. S hood & acrylic covers and Lower half of the machine complete with Mat finished S.S. covers. This reduces the noise level and to & fro flow of dust between compression zone and compression room.
- Antivibration mounting provided, Hence no foundation required.
- Large paint free compression zone with plenty of space for ease of cleaning & Accessibility of components.
- Electronic counter showing total output, output per minute. & RPM is provided on the control panel.
- Operating system same as standard machine and it is easy to operate and maintain the machine.
- Most of the components are common with the standard machine which reduces inventory of spares.
- Higher Turret speed upto 42 RPM gives higher output than standard machine.
- Upper punch dust covers are provided to aviod collection of dust on punch head.
- Maximum Turret RPM for GMP Machine is 42.

RIDDHI PHARMA MACHINERY LTD.

Manufacturer & Exporters of Pharmaceutical Machineries

OFFICE: MUMBAI Om-Shivam Bldg.,2, Tarun Bharat, Sahar Rd, Andheri(E), Mumbai-400 099.Maharashtra, India. Tel.: 91-22-2832 1332 / 2838 9432 Fax: 91-22-2838 2479 Res: 91-22-2683 3777, Mob: 98210 18579 Whatsapp: 919821018579 Skype: hiten.k.shah

E-mail: hitenkshah@hotmail.com/ hitenks@vsnl.com **Website:** www.riddhipharma.com **FACTORY:** Riddhi Estate, Gulab Nagar, Opp. AEC, Amraiwadi, Ahemdabad-380 026. India.



SINGLE ROTARY TABLET PRESS

Technical Specifications for Single Rotary Tabletting Machine																
Machine Model	Туре	Max Tab. Dia.	Max Depth Of FillW	Turret RPM		Output Per Min.	OP. Force In Tons	Tooling Type v		Motor HP	Case Dimensions L x W x H			Net Wt.	Gross Wt.	Remarks Max Length
				Min.	Max	Nos.*	1111016	Die	Punch	440V	Cms.			Kgs.	Kgs.	of DS.
RD D3-16	STD	23	20	16	30	260-480	10	D	D	2	120	127	183	980	1070	Max. Da. can be
RD D3-16	GMP	23	20	16	42	260-672	10	D	D	2	140	140	216	1070	1250	25.4 mm on demand
RD D3-20	STD	23	20	16	30	320-600	10	D	D	2	120	127	183	980	1070	25 mm
RD D3-20	GMP	23	20	16	42	320-840	10	D	D	2	140	140	216	1070	1250	for DS tablet
RD D3-23	STD	16	20	16	30	370-690	10	В	D	2	120	127	183	980	1070	19 mm for DS
RD D3-23	GMP	16	20	16	42	370-966	10	В	D	2	140	140	216	1070	1250	Tablet

^{*} Subject to the characteristic of the product

W 30 mm on demand at extra charges. v

D type .. punch/die dia 1" / 1 $^{1/2}$ $\rm v$

B Type Punch/die dia 3/4" 13/16"

OUR RANGE OF PRODUCTS

FOR COMPRESSION DEPT. : v 10 Stn. Mini Press GMP

- v 12 Stn. Slugging Machine STD/GMP v 16/20/23 Stn Single Rotary STD/GMP v 27/33/35/41/45 Stn Double Rotary STD/GMP
- v 29/35/37/45 Stn Press GMP
- v Dies & Punches
- v 150/300 CFM Dust Extractor Units
- v Dedusting & Deburring Unit.

FOR GRANULATION DEPT.:

- v Mass Mixer
- v Multi Mill
- v Communiting Mill
- v Tray Drier v Fluid Bed Drier
- v Roll Compactor

FOR COATING DEPT.:

- v 12/18/24/30/36/42/48/60/72 inch Coating machine
- v Polishing Pans for above M/c.
- $v \; Sprayin\bar{g} \; System$

OR LIQUID / OINTMENT DEPT.:

v Colloid Mill-plain and water iacketted

RIDDHI PHARMA MACHINERY LTD.

Manufacturer & Exporters of Pharmaceutical Machineries

OFFICE: MUMBAI Om-Shivam Bldg., 2, Tarun Bharat, Sahar Rd, Andheri(E), Mumbai-400 099. Maharashtra, India. Tel.: 91-22-2832 1332 / 2838 9432 Fax: 91-22-2838 2479 Res: 91-22-2683 3777, Mob: 98210 18579 Whatsapp: 919821018579 Skype: hiten.k.shah

E-mail: hitenkshah@hotmail.com/ hitenks@vsnl.com Website: www.riddhipharma.com FACTORY: Riddhi Estate, Gulab Nagar, Opp. AEC, Amraiwadi, Ahemdabad-380 026. India.