



Pill Tablet Press Machine TDP-5



Pill Tablet Press Machine is used for pressing different kinds of granular raw materials into round tablets. Pill Tablet Press Machine is applicable to trial manufacturing in lab or batch production in small amount different kinds of tablets, sugar pieces, calcium tablets and tablets of abnormal shape.

Pill Tablet Press Machine main application:

The machine is small desk-top tablet press, which can press various tablets continuously by means of manual operation.

With simple structure, convenient operation, Pill Tablet Press Machine has high degree of automation. Pill Tablet Press Machine can be used for small batch or laboratory production of a variety of medical tablets, sugar tablets, calcium tablets and so on.

The machine can press varieties of small particles and coarse powder into round or irregular tablets. It is not suitable to press large granular solids, wet powder or fine powder.

What is the TDP-5 Single punching machine suitable for suppressing:

TDP-5 single punch tablet press is suitable for pressing traditional Chinese medicine tablets, western medicine tablets, chewable tablets, calcium tablets, buccal tablets, spirulina tablets, milk tablets, candy tablets, effervescent tablets, medlar pills (health balls), aromatic tablets, taste blocks, disinfection tablets, catalysts, fertilizer tablets, pesticide tablets, button batteries, electronic components, metallurgical powder ceramics and so on.

Technical Parameter:

Model	TDP-5
Max. Pressure (kn)	50
Max. dia. of tablet (mm)	20
Max. Depth of fill (mm)	16
Thickness of largest tablet (mm)	6
Production capacity (pcs/h)	6000
Motor (kw)	0.75
Overall size(mm)	700*400*800
Net weight(kg)	125