



## ***RJWJ-200 Automatic Softgel Encapsulation Machine and Soft Gelatin Capsule Machine***



RJWJ-200 type soft capsule machine rolling mode is successful development on the basis of absorbing at home and abroad advance technology. This is roll mode soft capsule machine. The mold's length is 200 millimeters which have the new Row-hole approach to meet the needs of manufactures of pneumatic fixtures and micro-lubrication technology in China is a unique. Soft capsule machine can be oil. Suspension or paste into the quantitatively, produce different size, shape color of soft capsules. The soft capsules with fast and high bioavailability, elegant appearance, easy to swallow and difficult to oxidation, easy to storage and carrying it, etc. More and more consumers like it. Has been widely used in medicine, health products, cosmetics and other industries.

RJWJ-200 Series soft capsule machine is one of our company's new products which is automated production lines. The machine use new technology and materials, well-designed, precision manufactured. The molds uses Gas pressure to avoid disoperation of the damage caused by mold. Done to reduce the consumption.

Control box using the German Siemens and Siemens PLC programmable control true color touch screen, power control components are imposed to control the electronic control system, stability, performance, well, soft capsule machine in order to ensure the targets.

- 1). SIEMENS PLC control used (Compressive. Long life)
- 2). With fault diagnosis function
- 3). Simulation of machine images, vivid clear at a glance
- 4). Simple intuitive operation

## Parameter

<b>Model</b>	<b>RJWJ-200</b>
Roll mode speed	0-5rpm
Capacity	62400/h 8#olive,
Roller die size	103*200mm
Single-piston feed rate	0-2ml
Load difference	$\leq \pm 2\%$
Filling precision	$\pm 1\%$
Consumption of lubricating of rubber	~300ml/h
Power supply	380V 50HZ
Total power	9kw
Dimension	2050*890*2000mm
Weight	1400kg