CARMAPHARM

NJYF-300C Automatic Hard Capsule Liquid Filling and Sealing Machine



NJYF-300C automatic hard-capsule liquid filling machine and sealing linkage line equipment are the first in-line equipment based on our own multi-year capsule filling technology, self-developed design and production of hard capsule filling liquid and sealing. And the device can complete the filling of hard capsules (solution, suspension, micro emulsion or hot solution) and sealing the capsule cap joint to make the capsule. The contents are sealed in the process of packaging, transportation and use, which improves the stability and safety of the product. It can provide new modes of administration, improve the drug delivery problems of many drugs, and fill the domestic and foreign pharmaceutical and health care. The blank of hard capsule filling liquid and sealing technology in the product industry is also the gospel of many pharmaceutical companies and health care products factories for the development of some (highly active drugs, poorly soluble drugs) preparation products.

Features

- 1. The equipment adopts man-machine interface control system, which automatically lacks capsules and lacks material alarm.
- 2. The capsule feeding machine has a vacuum positioning and separating mechanism, so that the probability of the capsule is over 99%. The fiber is automatically detected, and the capsule is missing in the upper and lower mold holes.
- 3. High-precision ceramic pump and servo motor dosing system are adopted to control the difference of the loading amount to about ±2%, and the dosage is accurate.
- 4. The sealing mechanism adopts the servo motor lifting device, which can be used to jog and fine-tune the rubber wheel to achieve the best effect of the sealing layer.

- 5. The sealing liquid tank is designed as a sealing structure, and is equipped with a heating and heat preservation system to ensure the viscosity of the glue liquid and effectively extend the glue liquid period.
- 6. The finished dry curing mold adopts hanging type and free rotation structure, which effectively improves the disassembly and unloading cleaning time.

Technical Specification:

| Model | NJYF-300C |
|----------------------------|--|
| Production Capacity | 18,000 capsules / hour |
| Sealable capsule filling | Solution, suspension, microemulsion |
| Mold hole number | 6 holes |
| Applicable capsule | 00-4# gelatin, plant capsule |
| Environmental requirements | Room temperature around 20C, humidity 40% -55% |
| Voltage | Three-phase 380V/220V 50Hz |
| Total power | 9KW |
| Compressed air | 0.6Mpa |
| Total weight | 2300KG |
| Dimensions | 3800*900*1900mm (length × width × height) |