

DJL-8z Desktop Electronic Counter



The equipment counts and fills medicine automatically but bottles should be positioned on manual.

If the kind or filling quantity of medicine is changed, you only need to do simple setup on the touch screen.

This equipment is mainly used for counting and filling of tablets (including abnormal tablets), capsules, gels (including transparent gels and opaque gels) pills and similar solid granules. Using our own original high anti-dust sensor technology, the machine works accurately with good stability in heavy dust conditions.

Feeding evenly with quiver, special patent filling part protects materials from harm.

Technical Parameter:

Applicable	0-5# capsules, soft gels, ϕ 5.5-15 normal tablets, special shape tablets, coating tablets, ϕ 3-12 pills
Capacity	10~30 BPM
	4500~13500 bottles/shift
Working range	Adjustable from 2-9999 per bottle
Power supply	220V/50Hz 0.6kw
Size	660*1280*780mm
Weight	120kg
Remarks	1. Capacity is subject to the type and size of the bottle and tablet, as well as the quantity filled per bottle; 2. *BPM=Bottles per minute