

LHSS SERIES LABORATORY HIGH SHEAR MIXER



LHSS LABORATORY HIGH SHEAR MIXER, ONE OF THE COMMONLY USED EQUIPMENTS IN LABORATORY SOLID PREPARATION PRODUCTION, IS WIDELY USED IN PHARMACEUTICAL, CHEMICAL AND FOODSTUFF INDUSTRIES. THIS EQUIPMENT HAS FUNCTIONS OF MIXING AND GRANULATING.

FEATURES

- 1. FOLLOW THE REGULATIONS OF CGMP, PIC/S GMP, FDA
- 2. MOVABLE DESIGN
- 3. CHANGEABLE GRANULATING BOWL

ADVANTAGES

- 1. SMALL FOOTPRINT
- 2. R&D AND SCALE UP
- 3. GRANULATION CAN BE TIME, VOLUME AND POWER CONTROLLED

MAIN TECHNICAL PARAMETER

| Model | | LHSS5 | LHSS10 | LHSS15 |
|-------------------------|--------|-----------|-----------|-----------|
| Н | | 1100 | 1355 | 1400 |
| H1 | | 930 | 1020 | 1000 |
| H2 | | 755 | 755 | 700 |
| H3 | | 740 | 760 | 690 |
| Ĺ | | 1090 | 1165 | 1720 |
| L1 | | 920 | 920 | 1300 |
| G | | 550 | 550 | 650 |
| Volume | L | 5 | 10 | 15 |
| Feed capacity | kg | 1~2 | 2~4 | 3~6 |
| Mixing Power | kW | 1.1 | 1.1 | 2.2 |
| Chopper Speed | rpm | 90~900 | 90~900 | 50~290 |
| Granulation Power | kW | 0.75 | 0.75 | 1.1 |
| Impeller Speed | rpm | 1000~2800 | 1000~2800 | 1000~2000 |
| Comp-air Pressure | MPa | 0.6 | 0.6 | 0.6 |
| Comp-air Consumption | m³/min | 0.05 | 0.05 | 0.4 |
| Machine weight | t | 0.12 | 0.25 | 0.3 |