



Rotary Chicken Essence Tablet Press Machinery ZP-17B



The type of Rotary Chicken Essence Tablet Press Machinery is ZP-17B. The ZP-17B rotary tablet press is a rotary tablet press newly developed by our company to make granular materials into tablets of various shapes.

Rotary Chicken Essence Tablet Press Machinery maintenance

1. Regularly check the mechanical parts of the Rotary Chicken Essence Tablet Press Machinery, and carry out 1-2 times a month. Check whether the movable parts such as the worm gear, the worm, the bearing, the pressure wheel, the crankshaft, the upper and lower rails are flexible and wearable, and the defects should be timely. Repair use.
2. After using the Rotary Chicken Essence Tablet Press Machinery, remove the remaining powder and brush off the residual powder in the machine name. If the deactivation time is long, all the dies must be removed, and the machine should be wiped clean and the machine is smooth. Apply anti-rust oil and cover with a cloth.
3. The die maintenance of the Rotary Chicken Essence Tablet Press Machinery should be placed in a covered metal box, so that the die is completely immersed in the oil, and should be kept clean, not to cause rust and bruise. Customize the iron box as much as possible to fit each specification. Box, can avoid misuse and help to grasp the damage

Technical Parameter:

Dies (sets付)	17
Max. Pressure (kn)	120
Max. dia. of tablet (mm)	36
Max. Depth of fill (mm)	18-30
Max. Thickness of largest table (mm)	8-15
Turret speed (r/min)	10-25
Length of up and down punch pole (mm)	45
Length of up punch pole (mm)	175
Length of down punch pole (mm)	180
Diameter of middle mould (mm)	52
Production capacity (pcs/h)	11400-25500
Motor power (kW)	7.5
Overall size (mm)	900*800*1640
Weight (kg)	1500