



HM VL-CH series Vial Capping Machine



Usage: crimped seal vial with Alu-cap or Alu-plastic compound cap

Suitable: 2-50ml full stopper vial or other glass bottle;

Working principle: automatic feed vial and cap pull on cap, single-knife crimp cap, then count and discharge

Range: capping rate $\geq 99\%$; feeding rate $\geq 99\%$

Model: HM VL-CH400/HM VL-CH250

Capacity: 7ml vial, 400/250bpm; 6000-25000/5000-18000bph.