

SHH-KZ series Automatic Pharmaceutical Glass Injectable Vial Liquid Filling Sealing Washing and Capping Machine



Features

- 1. The compact line realizes single linkage, continuous operation from washing, sterilizing & drying, filling & stoppering and capping. The whole production process realizes cleaning operation; protects products from contamination, meets the GMP production standard.
- 2. Full servo control.
- 3. Transparent self-lifting protective cover with humid air outlet, electric screw control, safe and easy to maintenance.
- 4. For customers' liquid medicine and filling accuracy requirements, the ceramic pump filling system is selected, which can effectively ensure the filling accuracy and can be flexibly switched.
- 5. The stoppering form of inserting while rotating can effectively ensure the stoppering effect.
- 6. Capping machine: No Vial No Capping, No Stopper No Capping, Vacuum absorb aluminum scrap device.

Technical Parameters:

| Model | Production line | Suitable size | Output(max) | Power | Net weight | Overall size |
|---------------|------------------|---------------|-------------------|-------------|---------------|-----------------------|
| SHH-KZ I | CLQ 40 | 2-25 ml | 6000-12000 pcs/h | 69.8 KW | 7500 Kg | 9930×2500×2340 mm |
| | RSM 620/4250 | | | | | |
| | KGF 8 | | | | | |
| SHH-KZ II | CLQ 60 | 2-25 ml | 8000-18000 pcs/h | 85.8 KW | 8000 Kg | 10830×2500×2340 mm |
| | RSM 620/6000 | | | | | |
| | KGF10 | | | | | |
| SHH-KZ III | CLQ 80 | 2-25 ml | 10000-24000 pcs/h | 123.8 KW | 8100 Kg | 10830×2500×2340 mm |
| | RSM 900/10000 | | | | | |
| | KGF 12 | | | | | |