

SED-J Series High Speed Automatic SUS316 Hard Gelatin Capsule Filling Machine Pharmaceutical Equipment



This machine is designed further improved on the base of original full auto. Capsule filling machine from home and abroad, with high technology and exclusive performance, its function can reach international leading it is an ideal equipment for capsule and medicine in pharmaceutical industry.

Main Feature

- 1). Good-looking appearance, exquisite workmanship, ease of operation, simplicity of use.
- 2). Stowage seat and measuring plate are designed as one unit to make measuring plate and stowage rod without deviation, phenomenon avoid friction phenomenon between stowage rod and measuring plate, improve its precision highly, further more, it prolongs the machine's life greatly.
- 3). Ineligible capsule can be eliminated automatically (qualified rate not included), the medicine in the capsule can be recycled and be reused thus it can increase greatly economic efficiency.
- 4). Simplicity and convenience of dismantling, installation and clean, various models of mould can be replaced each other, the mould of 80 model and 1000 model as well as 1200 model can be replaced mutually on the same machine to meet different requirement of capacity.
- 5). Dust collector and vacuum pipe as well as waste air pipe are installed in the inner of the machine to avoid the air pipe becoming hard, bard broken and leakage etc phenomenon, it is more convenient to clean the platform. Furthermore, it accords with the requirement of GMP that medicine can not contact with organic material.

Product Parameters

Model	SED-200J	SED-400J	SED-600J
Max.output(capsules/min)	200	400	600
Filling dosage	Powder, Pellets, Tablet	Powder, Pellets, Tablet	Powder, Pellets, Tablet
Capsule size	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E
Total power	220V/380V 3.7Kw	220V/380V 3.7Kw	220V/380V 4.9Kw
Water supply	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa
Dust collector	160m3/h	160m3/h	160m3/h
Overall dimension (mm)	720*680*1700	720*680*1700	930*790*1930
Weight (kg)	700	700	800

Model	SED-800J	SED-1000J	SED-1200J
Max.output(capsules/min)	800	1000	1200
Filling dosage	Powder, Pellets, Tablet	Powder, Pellets, Tablet	Powder, Pellets, Tablet
Capsule size	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E
Total power	220V/380V 4.9Kw	220V/380V 5.75Kw	220V/380V 5.75Kw
Water supply	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa
Dust collector	160m3/h	160m3/h	160m3/h
Overall dimension (mm)	930*790*1930	1020*860*1970	1020*860*1970
Weight (kg)	800	900	900

Model	SED-1500J	SED-2000J	SED-2500J
Max.output(capsules/min)	1500	2000	2500
Filling dosage	Powder, Pellets, Tablet	Powder, Pellets, Tablet	Powder, Pellets, Tablet
Capsule size	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E
Total power	220V/380V 6.45Kw	220V/380V 6.45Kw	220V/380V 6.45Kw
Water supply	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa
Dust collector	160m3/h	160m3/h	160m3/h
Overall dimension (mm)	1020*1050*1970	1020*1050*1970	1150*1050*1970
Weight (kg)	900	1500	1500

Model	SED-3200J	SED-3800J	SED-7500J
Max.output(capsules/min)	3200	3800	7500
Filling dosage	Powder, Pellets, Tablet	Powder, Pellets, Tablet	Powder, Pellets, Tablet
Capsule size	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E	Capsule 00# -5#/Safety capsule A-E
Total power	220V/380V 9.55Kw	220V/380V 9.55Kw	220V/380V 13.5Kw
Water supply	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa	500L/h,0.2-0.10Mpa
Dust collector	160m3/h	160m3/h	160m3/h
Overall dimension (mm)	1150*1280*2060	1150*1280*2060	1700*1280*2150
Weight (kg)	2400	2400	3800