

## Counting DJD 1A



The single-file equipment counts and fills capsules, soft gels, normal tablets, special shape tablets, coating tablets, pills and most solid granules in small bags. Counts by infrared and sensor that doesn't wear or break the product and is compatible with all kinds of materials. Can accomplish bag forming, counting, filling, bag sealing and printing at the same time.

Modularized design for very easy cleaning and disassembling, only need to change a couple of parts when shifting between different products.

### Application:

Auto bag forming, counting, filling and bag sealing and printing at the same time.

Suitable for filling in small bags with back sealing, 3 sides sealing, 4 sides sealing (according to user's requirements).

### Technical Parameter:

|                            |                                      |
|----------------------------|--------------------------------------|
| Length of bag (mm)         | 50-150                               |
| Width of bag (mm)          | 50-100                               |
| Capacity (bags per minute) | 20-45                                |
| Product Sensor (Infrared)  | To count and detect breakage product |

|                          |  |
|--------------------------|--|
| Working Range            | lebih dari atau sama dengan 2 (adjustable)       |
| Counting system          | Nantong Hengli                                   |
| Speed Control            | T.W.T  |
| Power Supply             | 220 v ; 50 Hz                                    |
| Power Consumption        | 1.5  |
| Weight                   | 300  |
| Dimension (L x W x H) mm | 650x1080x1630                                    |
| Main Motor               | GPG (China)                                      |
| HMI                      | HITECH (Taiwan)                                  |
| Electrical Part          | Schneider (France)                               |
| Emergency stop           | Stop main machine in case of emergency situation |