

## *HM VL-RW series Rotary Ultrasonic Vial Washing Machine*



Usage: rotary movement, Ultrasonic washing, water and compressed air multi-stages wash/clean and dry.

Vial breakage:  $\leq 0.1\%$

Clarified vial:  $\geq 99\%$

Model: HM VL-RW400/300/150

Washing heads: 80/60/40 pieces

Capacity: 6000-24000 vials/h, 400 vials/min for 7ml vial;  
5000-18000 vials/h, 300 vials/min for 7ml vial;

3000-9000 vials/h, 150 vials/min for 7ml vial.