



— Pharma Process & Packaging

www.saintyco.com

Uphill Tablet Deduster

C&C200F



CATALOGUE

Uphill Tablet Deduster

C&C200F Uphill Tablet Deduster is a special equipment to wipe out burr and dust of tablet whilst going up way. It can be connected with tablet press and metal detection systems. 316L stainless steel are used for all metal parts contacting tablets. Design of GMP.



Features:

- Capacity: 800,000 tablets/hr. (ref. Ø6mm round tablet)
- Body and cover of the machine take building-block design. Height of outlet is easy to adjust;
- 360° circle adjustment of inlet and outlet to make production more convenient and save factory room;
- Fast-joint system to ensure safe and quick disassembly and assembly;
- Outstanding efficiency can make mass production of $\Phi 3$ - $\Phi 35$ mm all-kind-shape tablet;
- Parts in contact with medicine all made of 316L tainless steel. No blind spot during whole machine cleaning. Meet cGMP demands;
- New design motor driver to make machine run more smoothly with lower noise.



C&C200F+Metal Detector

Technical Specifications	C&C200F
Tablet Diameter (mm)	Ø3-Ø25
Max Output Capacity (t/h)	800000(Ø6)
Inlet height (mm)	790-970
Outlet height (mm)	1110-1290
Dedusting Distance (m)	7.5
Power Supply	1PH 220V/110V 50Hz/60Hz 0.1kw
Compression Air	0.1 m ³ /min 0.05MPa
Vacuum	2.7 m ³ /min -0.1MPa
Dimension (mm)	515 x 545 x(1345-1525)
Weight (kg)	80

SaintyCo also offer:

Encapsulation Solutions

Tablet Compression Solutions

Coating Systems

Cartoning Systems

Blister Packing Solutions

