



## *Pharma Vial Stoppering Filling Capping Labeling Machine*



Terminal sterilized and non terminal sterilized Water injection, lyophilized injection powder, IV solution and other medicine.

### **Main characteristic:**

1. Achieve full time weighing detect of Whole process and automatically feedback the information to adjust filling volume. Working principle of domestic existing type is: firstly take out bottles, weighing then filling, then get bottle for weighing again then take back for closing, this may cause cross-contamination and omission. But our on-line whole process of weighing detection can avoid such kind of risk.
2. Using straight-line stoppering, will not disrupt laminar flow, can be fully avoided rotation filling block of laminar flow and mixed flow caused by rotation, turbulence, avoid cross contamination risks, so meet the new GMP requirements.
3. "0 risk" designing according to the most advance international new concept, there is no mechanical transmission chain on machine table, no dead area, and machine table completed separated from transmission box, so that no friction dusts and oil pollution, which is the pioneer in China, and it's the upgrade product of the same products.

**Technical Parameters:**

Product Name	Vial Filling-Sealing Machine
Applicable vial	2-100ml
Production capacity	2ml:150 Vial/min~400 Vial/min
Filling head	4~12
Filling error	±1%
Broken rate	≤0.1%
Stoppering rate	≥99.9%
Dimension	2,4*1,7*2,6 (m)
Weight	2700 Kg~3000 Kg