

## HSD Series Laboratory Mixer



It is mainly used in laboratory solid preparation for mixing of granules and granules, granules and powder, powder and powder. Meanwhile, it is widely used in pharmaceutical, chemical and foodstuff industries.

### Features

- Movable design
- Detachable bin design for various batch size application
- Efficient blending

### Advantages

1. Suit for R&D
2. High efficiency mixing
3. Small footprint

### Parameters

Model		HSD5	HSD10	HSD15	HSD20	HSD25
Blending speed	rpm	3~25				
Max.bin Volume	L	5	10	15	20	25
Max.loading volume	L	4	8	12	16	20
Max. loading weight	kg	2.5	5	7.5	10	12.5
Total power	kw	0.37	0.37	0.37	0.37	0.55
Total Weight	kg	150	160	170	180	185

Model		HSD30	HSD50	HSD75	HSD100
Blending speed	rpm	3~25			3~20
Max.bin Volume	L	30	50	75	100
Max.loading volume	L	24	40	60	80
Max. loading weight	kg	15	25	37.5	50
Total power	kw	0.55	0.55	0.75	0.75
Total Weight	kg	190	200	230	260