



Pharmaceutical Rotary Tablet Press Machine ZP-5/7/9B



Features:

1. The peripheral nozzle of the Pharmaceutical Rotary Tablet Press Machine is fully enclosed, the material is made of stainless steel, and the internal countertop is made of stainless steel, which can maintain the surface gloss and prevent cross-contamination, in line with GMP.
2. With the transparent protective cover, you can clearly see the operation of the Pharmaceutical Rotary Tablet Press Machine. It can be fully opened for easy cleaning and maintenance.
3. The speed of the Pharmaceutical Rotary Tablet Press Machine is stepless speed regulation, which can be adjusted freely.
4. Realize electromechanical integration, the tablet press adopts touch mode key and display screen (speed, output, fault), etc.
5. The transmission system is sealed in the oil tank below the main body of the machine. It is a separate component that is safely separated. It does not pollute each other, and the transmission shaft is infiltrated into the oil pool, which is fully lubricated, and the noise and wear are reduced.
6. The Pharmaceutical Rotary Tablet Press Machine is equipped with a vacuum cleaner, which absorbs the dust from the indoor air and absorbs it.

Technical Parameter:

Type	ZP-5B	ZP-7B	ZP-9B
Dies	5	7	9
Max.pressure(kn)	60	60	60
Max.Dia.of Tablet(mm)	20	20	20
Max.Depth of Fill(mm)	15	15	15
Max.Thickness of Tablet(mm)	6	6	6
Turret Speed(r/min)	30	30	30
Production Capacity(pc/h)	9000	12600	16200
Motor(kw)	2.2	2.2	2.2
Overall Size(mm)	480*630*1100	480*630*1100	480*630*1100
Net Weight(kg)	360	360	360