



SLAG-2A Pharmaceutical Ampoule Filling Sealing Machine



Introduction

1. Filling and sealing machine for ampoule Model SLAG-2A is equipment specifilized in filling and sealing 1ml to 2ml, 5mlto10ml, and 15-20ml bottle in pharmaceutical works.
2. This machine is equipped with automatic stop-filling device which can make filling case when there is no ampoule into filling station.
3. Suitable for pharmaceutical ampoule filling and sealing drawing, the local station for the 2 bottles for one stepper drive system, each with 2 bottles as a group, were sent to the various processes point to complete the delivery bottles, before nitrogen, drench, after nitrogen, preheated, drawing sealing processes.
4. Using hydrogen generator, the interface is fully matched. Drive rack stainless steel, pressure bottle tied tight fit.

Specification:

Filling Nozzles	Applicable Volume	Yield Capacity(p/h)
2 Nozzles	1-2ml	2500-3500
	5-10ml	2000-3000
	15-20ml	1200-1500
4 Nozzles	1-2ml	7000-9000
	5-10ml	6000-8000
6 Nozzles	5-10ml	3000-5000
8 Nozzles	1-2ml	1600-1800