



McBRADY ENGINEERING, INC.



<http://www.mcbradyengineering.com/>

1251 S. LARKIN AVE.
JOLIET, IL 60436



EXTERIOR VIAL WASHER/DRYER

The Exterior Vial Washer/Dryer will wash, rinse and dry the exterior of filled and capped vials/ampoules. The entire exterior of each vial is exposed for cleaning as it passes through our patented “Zig-Zag” guide pattern. A low-pressure wash and rinse are used to keep the cap area dry. The option of multiple passes (1 to 10) including choice of rinse sequence makes the Exterior Vial Washer/Dryer the most versatile in the industry. Capable of reaching speeds of 150vpm - 500vpm (depending on vial size), this machine is the perfect addition to any pharmaceutical line.



Special Features

- Multiple rinses/dry passes allow for **“bone dry” vials**
- **No Change Parts.** Changeovers for the single pass washer are **quick and easy**, simply adjust the rails for new vials
- **No cap wetting.** The low pressure wash gently cleans the exterior of vials, while leaving the caps spotless
- **High Speed Cleaning.** The 10-pass vial washer can reach speeds of up to 500vpm (depending on vial size)

MEMBER



SPECIFICATIONS

SINGLE PASS DIMENSIONS

A: Overall Width
24"

B: Overall Length
6' → 10'



MULTIPLE PASS DIMENSIONS

Pass Qty.	Width (A)	Length (B)	Speed
4 Pass	24"	48"	150
6 Pass	32"	48"	200
8 Pass	43"	58"	350
10 Pass	54"	68"	500

**Speed in vials per minute*

**Actual speed determined by vial size*

