

LD Series Mixing Bin



LD Series Mixing Bin could be used in HLD Series Bin Mixer and SLD Series Bin Lifter, all using the same bin in different technology locations. Thus avoiding the frequent material transfer dust generated and cross-contamination, so that drug production in accord with GMP requirements.

Use high quality austenitic stainless steel SUS304 or SUS316L, smooth inner surface (Ra<0.4 μ m) no dead ends, easy to clean. Bin hoper lid with silicone rubber seal, discharge valve with antimisuse means reliable sealing. Polyurethane casters, leaving no trace, no damage to the ground floor.

Technical Parameters:

| Model | LD400 | LD600 | LD800 | LD1000 | LD1200 | LD1500 | LD1800 | LD2000 |
|--------|---------|---------|---------|---------|---------|---------|---------|---------|
| Volume | 400 L | 600 L | 800 L | 1000 L | 1200 L | 1500 L | 1800 L | 2000 L |
| Α | 800 mm | 900 mm | 1000 mm | 1000 mm | 1200 mm | 1200 mm | 1200 mm | 1200 mm |
| Н | 1150 mm | 1400 mm | 1500 mm | 1553 mm | 1640 mm | 1835 mm | 1990 mm | 2230 mm |
| H' | 120 mm |