



## *High Quality Mini Single Punch Tablet Press DP-12/25*



### **High Quality Mini Single Punch Tablet Press main feature:**

The fully closed machine which High Quality Single Punch Tablet Press avoids dust cross contamination and conforms to GMP requirements.

The safety protective device adopts more reasonable user-friendly design.

Single Punch Tablet Press has wider operation space. With prominent bearing capacity, Single Punch Tablet Press For Pharmaceutical Machinery is suitable for production of a wider range of tablets.

### **High Quality Mini Single Punch Tablet Press Main purpose:**

Researched and developed successfully by absorbing advanced technologies from abroad, the High Quality Mini Single Punch Tablet Press is a continuous automatic tablet press used for pressing powder and granular raw materials into tablets in such industries as medicine, chemical industry, electronics and food.

With high utilization rate of materials, reasonable structure, convenient operation, low volume, high pressure, low noise, little difference in tablet weight, continuous adjustment of filling depth of materials and tablet thickness, all performance indexes of this machine reach the level of an imported one. Fully conforming to requirements for pharmaceutical sanitation in China, Mini Single Punch Tablet Press is a preferred machine for scientific research institutes, pharmaceutical plant process rooms and hospital preparation rooms.

Mini Single Punch Tablet Press adopts many novel designs during its designing process, so as to ensure that the production process is safer and more standard.

**Technical Parameter:**

<b>Model</b>	<b>DP-12</b>	<b>DP-25</b>
Punch quantity(set)	1	1
Max.tablet diameter(mm)	12	25
Max.output(tablets/minute)	60	60
Max.main pressure(kn)	50	50
Max.filling depth(mm)	20	20
Max.tablet thickness(mm)	8	8
Motor power(kw)	220V/50Hz,1.1KW	220V/50Hz,1.5KW
Overall size(mm)	550*450*750	580*500*830
Machine weight(kg)	100	150