

台中市北區大仁街 81 號

81 Tajen St, N, Dist., Taichung, Taiwan

Tel:886-4-23227637 Email:shingtse@ms23.hinet.net Http://www.shingtse.com

# HIGH SPEED ROTARY TABLET MACHINE TCS-HRT31E-SS TCS-HRT39E-SS

(ONE-LAYER TABLETS)

#### **FUNCTIONS**

- Machine is ideal for making round and irregular shape tablets with high hardness, high precision and high production.
- Applicable to food, medicine, chemical, magnets & electronic components.

#### **FEATURES**

- Electrical control digital display meter adjusted device is for adjustment in hardness, thickness and weight of pastilles.
- With dual pressing wheel system: machine installs pre-pressing wheels and main-pressing wheels.
- The structures of two sets of pressing wheels are integrated casting to enhance the precision of pastilles.
- Tableting area separates from rotary gear to reduce noise (it is a split-type).
- Large hopper and new three-phase feeding design for easy assembly and disassembly.
- Can use both optional fixed amount feeding device and forced feeding device to easily assemble and disassemble.
- Easy operation & cleaning as well as meet PIC/S GMP standards.





台中市北區大仁街 81 號

81 Tajen St, N, Dist., Taichung, Taiwan

Tel:886-4-23227637 Email:shingtse@ms23.hinet.net Http://www.shingtse.com

# HIGH SPEED ROTARY TABLET MACHINE TCS-HRT31E-SS TCS-HRT39E-SS

(ONE-LAYER TABLETS)

#### **SPECIFICATION**

| SPEC/MODEL                              | HRT31E-SS                | HRT39E-SS         |
|---|--------------------------|-------------------|
| TABLET TYPE                             | ONE LAYER                |                   |
| NUMBER OF STATION                       | 31                       | 39                |
| PUNCH TYPE                              | D                        | D                 |
| MAX. DIA. OF TABLET                     | 25 mm                    | 16mm,19mm(caplet) |
| DEPTH OF FILLING                        | 2 ~ 19 mm                |                   |
| PRODUCTION CAPACITY(TABLET/HR)          | 22,300 ~ 89,200          | 28,000 ~ 112,000  |
| RPM OF TURRET                           | 12 ~ 48 RPM              |                   |
| MAX. MAIN PRESSURE                      | 14,000 KGS               |                   |
| MAX. PRE-PRESSURE                       | 14,000 KGS               |                   |
| OVERALL LENGHT OF UPPER & LOWER PUNCHES | 133.6 mm                 |                   |
| DIAMETER OF UPPER & LOWER PUNCHES       | 25.4 mm                  | 25.4mm            |
| THICKNNES OF DIE                        | 23.81mm                  | 23.81mm,          |
| DIAMETER OF DIE                         | 38.1 mm                  | 30.16 mm          |
| POWER                                   | 3 HP 6 pole, 5 HP 4 pole |                   |
| MACHINE SIZE (L x W x H)                | 1,160 × 1,060 × 1,980 mm |                   |
| PACKING SIZE (L x W x H)                | 1,350 × 1,300 × 2,300 mm |                   |
| MACHINE TOTAL WEIGHT                    | 2,700 KG                 |                   |

Remarks: Output is available according to powder character and the evenness of powder granules.

The above figures are for reference only; the specifications can be modified by manufacturer as needed.

<sup>\*</sup>Punches and dies are available according to purchasers' requirements.



台中市北區大仁街 81 號

81 Tajen St, N, Dist., Taichung, Taiwan

Tel:886-4-23227637 Email:shingtse@ms23.hinet.net Http://www.shingtse.com

### **HIGH SPEED ROTARY TABLET MACHINE**

TCS-HRT31E-SS TCS-HRT39E-SS

#### **Feature explanation**

#### Pre-& Main-pressures and weight adjustment

- Use humanity electronic digital meter to display precisely, push button & labor-saving adjustment design.
- It is easy to control tablet production data and record QC efficiently.
- Can input the recorded datum to save the tableting time when the sizes of powder granule are the same and under the same working condition.
- With the brake motor transmission and mechanism regulated device to increase the stability of hardness
  of tablets efficiently.

#### Integrated double pressing system

- Double pressing wheels: two sets of pre-pressing wheels and main-pressing wheels, it is applicable to special powders for producing pastille. (see below the picture)
- Two sets of pressing wheel structures: integrated casting (non-pillar type) to reduce the noise when tablets are made and also enhance the precision of tablets. (see below the picture)

#### **Split-type transmission mechanism**

- Tableting area separates from rotary gear (it is a split-type); spindle drives the turret directly that can reduce the noise.
- The spindle worm gear inside of gear oil box meshes with the copper alloy warm gear to keep the noise down and extend the machine's life.

#### Safety device of punches & dies

- The turret has hard chrome plated surface treatment.
- Each turret aperture has the same diameter size; it is applicable to round type and special tablet forming punches.
- Upper punches install an oil-proof loop and lower punches install an oil seal to prevent the tablets from the oil sludge pollution.
- Upper punches equip with the safety device to stop machine automatically when it is stuck; lower punches have the protective cover with a safety sensor when changing the punches.

#### **Track device**

- All tracks, which connect with upper and lower punches, are made of special steel.
- Installation of the auto oil filling device that the lower track can fill oil and lubricate automatically for the smooth work of lower punches.

#### Feeding design

- With a large hopper and new three-phase-split type input design for easy assembly and disassembly.
- Can use an optional forced feeding device which is driven by variable speed motor for easy assembly and disassembly.





台中市北區大仁街 81 號

81 Tajen St, N, Dist., Taichung, Taiwan

Tel:886-4-23227637 Email:shingtse@ms23.hinet.net Http://www.shingtse.com

### **HIGH SPEED ROTARY TABLET MACHINE**

TCS-HRT31E-SS TCS-HRT39E-SS

### **Feature explanation**

### **Easy operation & cleaning**

- Machine exterior and interior tableting areas are clad by stainless steel#304 or have the hard chrome treatment.
- User-friendly electrical control panel.
- Equip with the tablet counter, machine working timer and electrical appliance safety device.
- All of panels, tracks and punches are easy to clean, assemble and disassemble.
- Meet PIC/S GMP standards.



台中市北區大仁街 81 號

81 Tajen St, N, Dist., Taichung, Taiwan

Tel:886-4-23227637 Email:shingtse@ms23.hinet.net Http://www.shingtse.com

### HIGH SPEED ROTARY TABLET MACHINE

TCS-HRT31E-SS TCS-HRT39E-SS

**Optional Model explanation** 

#### Handwheel adjusting machines

• Tablet hardness, thickness and weight adjustment: using handwheel adjustment device (digital meter display).



### Forced feeding & fixed amount feeding machines

- Applicable to granules that have a bad fluidity or powders which need to be tableted.
- Driven by variable speed motor for easy assembly and disassembly.
- See the right picture.

#### **Touch screen & PLC machines**

• With cantilever and electrical control box panels.

#### **Options**

- Air pin breaker
- Powder forced filling device
- Powder forced filling device + powder fixed amount offering device
- Touch screen + PLC device
- Deduster TCS-I-SS
- Various punches & dies specifications and sizes can be customized

