

## DJL-16 Electronic Counter



This equipment is mainly used for bottle packing of tablets (including abnormal tablets), capsules, gels (including transparent gels and opaque gels) pills and similar solid granules.

Highly accurate counting.

Using our own original high anti-dust sensor technology, the machine works accurately with good stability in heavy dust conditions.

Feeding evenly with quiver, special patent filling part protects materials from harm.

No bottle under the filling hopper, no material in the supply hopper, or the photoelectric sensor box is abnormal, the machine makes alarm.

The special filling hopper avoids the materials jam, small mouth bottle can quickly bottling.

(Available as an option)

You can choose the advanced device of getting rid of powder and the device of getting rid of static electricity. (Available as an option)

Sound and light alarm. (Available as an option)

Detection options include: lack of products in the supply hopper, damaged tablet, weighing systems online. (Available as an option)

Fully compliant to GMP regulations.

**Technical Parameter:**

Applicable	00-5#capsules,soft gels $\phi$ 5.5-22normal tablets, special shape tablets, coating tablets $\phi$ 5.5-20pills and other solid granules
Capacity	20~80 BPM
	9000~36000 bottles/shift
Working range	Adjustable from 2-9999 per bottle
Power supply	220V/50Hz 1.2kw
Air supply	0.6Mpa, 36-45l/min
Size	1800*1600*1560mm
Weight	410 kg
Remarks	1. Capacity is subject to the type and size of the bottle and tablet, as well as the quantity filled per bottle; 2. *BPM=Bottles per minute