

HCX rotary-type ultrasonic bottle-washing machine



Purpose and brief

Local low-dose bottles are dedicated cleaning equipment, bottles by the conveyor belt Net immersion in water for ultrasonic cleaning rough, and then one by one by the screw into the rising block, the bottles and then were trapped out of the water bottle clamp to live, after flip after , respectively, by the back wash water rough - blowing compressed air purification - pure water or distilled water to wash sperm - blowing compressed air purification, complete bottle-washing job. Finally the bottle and then turned round bottles from the allocation set aside.

Main technique parameters

Production capacity: 250 ~ 400 bottles / min

Applicable specification:

1ml ~ 20ml Xilin bottle, ampoule bottles, oral liquid bottles

Circulating water flow: 0.6m³ / h

Purified Water consumption: 0.4 ~ 0.5m³ / h Pressure: 0.2 ~ 0.3MPa

Purification of compressed air consumption:

40m / h pressure: 0.2 ~ 0.3MPa

Ultrasonic Power: 1Kw

Power: 380V 50Hz three-phase four-wire system

Machine power: 5Kw