

DWJ-2000D-2T Seamless Soft Capsule Machine



DWJ-2000D-2T dripping pill production line, it's with 200-hole automatic dripping pill machine, pill collecting machine and rotary type screening machine. It adopts the way that combines 10 pcs or more dripping system (every dripping system is with 100 holes, it can also be separately controlled) to make quantity of the dripping heads more than 1000pcs, since intensive control is adopted in cooling, PLC, touch screen, materials preparation tank and other parts, the automatic degree is higher, compared with other many production lines with the same production capacity, the production cost and price will be reduced. It is with small occupied area, small investment, easy to control, high rate of finished products, to control more production lines with the labor for one machine, and many other outstanding advantages.

Technical Parameters:

working pressure	380V
power	16kw
configuration	200 holes, or 40-80 holes
dripping container volume	1500mlx2
liquid tank volume	60L
external air pressure	$\geq 0.5\text{Mpa}$
external dimension	2600x800x1980mm
weight	660kgs