

## ***SED-10PSL Automatic Ultra wide Channel Capsule Counting Machine***



SED-10PSL Automatic Electronic counting machine is special for counting tablets, soft gelatin, hard capsule and chewing gums, etc. This machine adopts 10 sets imported counting sensor for counting with 10 extra widen guiding way and one head. The advantage of the machine is no need change mold when you change the counting object, only adjusting the counting table height by Easy Adjust wheel. The machine comes with Touch screen for easy operation.

### **Feature**

- 1) Compatibility, high accuracy, wide application and use of the Company's original anti-dust photoelectric detection of high technology, can count Tablets (including special-shaped tablets), capsules, soft capsules (transparent, opaque), pill count to bottling;
- 2) Vibratory cutting drugs do not damage;
- 3) A high degree of intelligence, with no bottle is not the number of self-fault detection alarm and many other control functions;
- 4) Compact structure, stability and durability, and no regular replacement parts, small footprint, the use of low-cost;
- 5) To facilitate the replacement of species, can automatically take-off and landing, demolition cleared without any tools, convenient and quick;
- 6) Chinese dialogue control panel, can store 10 ~ 20 group of species parameters, adjusted to facilitate the use of special;
- 7) Fully automatic operation, with the Upper and Lower Union Road, intelligent control of processes, coordination, and no special operation, labor-saving devices.
- 8) Low prices, especially for small-volume, variety and more manufacturers use. Control chip imported to the United States, photoelectric sensors for Taiwan's imports.

## Detailed Parameters

| Model             | SED-10PSL                                                                                                                                                             |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Productivity      | 20-50bottles/min                                                                                                                                                      |
| Applicable target | Filled capsule(000-5#)<br>Soft Gelatin capsule<br>chewing gum<br>tablet<br>other solid object under machine range                                                     |
| Hopper Volume     | 20-50bottles(30-100 pcs/bottle)/min                                                                                                                                   |
| Filling Quantity  | 1-9999pcs/bottle (adjustable)                                                                                                                                         |
| Bottle Size       | Diameter: 30-100mm<br>Height:40-180mm                                                                                                                                 |
| Dimension(L*W*H)  | 1360*1530*1700mm                                                                                                                                                      |
| Weight            | 320 kg                                                                                                                                                                |
| Counting Size     | Pills with diameter of 37 or less are compatible with 000-5 capsules soft capsules with diameters of 3.5-37 plain tablets special-shaped tablets sugar-coated tablets |
| Compressed air    | 0.6Mpa 3.5L/min                                                                                                                                                       |
| Power             | 220V/50HZ 0.6KW                                                                                                                                                       |